

Invacare® LiNX

Intuitive Dual Control (IDC)

en **Remote
User Manual**

This manual **MUST** be given to the user of the product.
BEFORE using this product, this manual **MUST** be read and saved for
future reference.



Yes, you can.®

Contents

1 General	1
1.1 About This Manual	1
1.2 Warranty	1
1.3 Service Life	1
1.4 Limitation of Liability	1
1.5 Symbols in This Manual	2
1.6 General Safety Notes	2
1.7 Intended Use	4
1.8 Components	4
1.8.1 Labels on Product	5
2 Setup	6
2.1 Installing / Removing the IDC	6
2.2 Connecting IDC to Controls System	7
3 Usage	8
3.1 Power up IDC	8
3.1.1 Switching Control Between Attendant and Occupant	8
3.2 Controlling the Power Wheelchair with IDC	8
3.2.1 How to Steer the Power Wheelchair	8
3.3 “Emergency Stop” Function	9
4 Troubleshooting	10
4.1 Fault Diagnosis	10
4.2 Fault Codes and Diagnosis Codes	10
5 Maintenance	13

© 2024 Invacare Corporation

All rights reserved. Republication, duplication or modification in whole or in part is prohibited without prior written permission from Invacare. Trademarks are identified by ™ and ®. All trademarks are owned by or licensed to Invacare Corporation or its subsidiaries unless otherwise noted.

1 General

1.1 About This Manual

This document is a supplement to the product's user documentation.

This component itself does not bear a CE and a UKCA mark but is part of a product that complies with the Medical Device Regulation 2017/745, Class I and Part II UK MDR 2002 (as amended) Class I concerning medical devices. It is therefore covered by the product's CE and UKCA marking. See the product's user documentation for more information.

Only use this component if you have read and understood this manual. Seek additional advice from a healthcare professional who is familiar with your medical condition and clarify any questions regarding the correct use and necessary adjustment with the healthcare professional.

Note that there may be sections in this document, which are not relevant to your component, since this document applies to all available models (on the date of printing). If not otherwise stated, each section in this document refers to all models of the component.

Invacare reserves the right to alter component specifications without further notice.

Before reading this document, make sure you have the latest version. You find the latest version as a PDF on the Invacare website.

Previous product versions may not be described in this Manual's current revision. If you require assistance, please contact Invacare.

If you find that the font size in the printed document is difficult to read, you can download the PDF version from the website. The PDF can then be scaled on screen to a font size that is more comfortable for you.

For more information about the component, for example safety notices and recalls, contact your Invacare representative. See addresses at the end of this document.

In case of a serious incident with the component, you should inform the manufacturer and the competent authority in your country.

1.2 Warranty

The terms and conditions of the warranty are part of the general terms and conditions particular to the individual countries in which this product is sold.

1.3 Service Life

We estimate a service life of five years for this product, provided it is used in strict accordance with the intended use as set out in this document and all maintenance and service requirements are met. The estimated service life can be exceeded if the product is carefully used and properly maintained, and provided technical and scientific advances do not result in technical limitations. The service life can also be considerably reduced by extreme or incorrect usage. The fact that we estimate a service life for this product does not constitute an additional warranty.

1.4 Limitation of Liability

Invacare accepts no liability for damage arising from:

- Non-compliance with the user manual
- Incorrect use
- Natural wear and tear
- Incorrect assembly or set-up by the purchaser or a third party
- Technical modifications and / or unauthorised modifications
- Use of unsuitable spare parts

1.5 Symbols in This Manual

Symbols and signal words are used in this manual and apply to hazards or unsafe practices which could result in personal injury or property damage. See the information below for definitions of the signal words.



WARNING!

Indicates a hazardous situation that could result in serious injury or death if it is not avoided.



CAUTION!

Indicates a hazardous situation that could result in minor or slight injury if it is not avoided.



NOTICE!

Indicates a hazardous situation that could result in damage to property if it is not avoided.



Tips and Recommendations

Gives useful tips, recommendations, and information for efficient, trouble-free use.

Other Symbols

(Not applicable for all manuals)



UK Responsible Person

Indicates if a product is not manufactured in the UK.



Triman

Indicates recycling and sorting rules (only relevant for France).

1.6 General Safety Notes



WARNING!

Risk of Injury or Damage to the Power Wheelchair

Do not install, maintain or operate this equipment before you have read and understood all the instructions and all the manuals for this product and all other products that you use or install together with this product.

- Follow the instructions in the user manuals.



WARNING!

Risk of Serious Injury or Damage to the Power Wheelchair or Surrounding Property

Wrong settings can make the power wheelchair uncontrollable or unstable. An uncontrolled or unstable power wheelchair can cause an unsafe situation such as a crash.

- Performance adjustments must only be made by qualified technicians or by persons who completely understand the programming parameters, the adjustment process, the configuration of the power wheelchair and the capabilities of the driver.
- Performance adjustments must only be made in dry conditions.

**WARNING!****Risk of Injury or Damage to the Power Wheelchair**

Risk of unintended movement of the power wheelchair or seating system when loose personal belongings (e.g. jewellery, scarfs) become entangled around the joystick.

- Make sure that any loose items are clear of the joystick when your power wheelchair is powered up.
- Power off your power wheelchair immediately to stop any movement.

**CAUTION!****Risk of Injury or Damage**

Malfunctioning or unintended movement may occur if the IDC was exposed to forceful impacts such as wall collision or fall damage, or if it is visibly damaged. Internal parts might be broken.

- If unintended movement occurs, stop using the IDC and power wheelchair immediately and contact a qualified technician.

**WARNING!****Risk of Injury or Damage due to Electrical Shorts**

Connector pins on cables connected to the power module can still be live even when the system is off.

- Cables with live pins should be connected, restrained or covered (with non-conductive materials) so that they are not exposed to human contact or materials that could cause electrical shorts.
- When cables with live pins have to be disconnected, for example, when removing the bus cable from the remote for safety reasons, make sure to restrain or cover the pins (with non-conductive materials).

**CAUTION!****Risk of Injury due to Unintended Movement**

It is recommended that the power wheelchair, fitted with a Gyro module, has a drive function with disabled Gyro. If the power wheelchair is used in a moving vehicle (e.g. boat, bus or train) maybe the Gyro function is impaired and drive demands can result in unintended movement.

- When driving on a moving vehicle choose a drive function with disabled Gyro.
- If the power wheelchair does not have a drive function with disabled Gyro, contact your Invacare provider.

**CAUTION!****Risk of Injury from Hot Surfaces**

Remote module can get hot when exposed to strong sunlight for long periods.

- Do not leave power wheelchair in direct sunlight for long periods.

! NOTICE!

If you touch the connector pins, they can become dirty or they can be damaged by electrostatic discharge.
 — Do not touch the connector pins.

! NOTICE!


There are no user-serviceable parts inside any case.
 — Do not open or disassemble any case.

1.7 Intended Use


The IDC can be used to assign control of the power wheelchair from the occupant to an attendant. The IDC is not suitable for the wheelchair occupant himself. The attendant walks behind or to the side of the wheelchair, and can intuitively determine speed and direction of movement by pressing a handle with minimum force. Adjustments can also be made to the wheelchair using dedicated programming. Once the attendant activates the IDC, he takes over complete responsibility for the wheelchair and its occupant.

1.8 Components


Top Side



Bottom Side



1.8.1 Labels on Product



	(A)	Position of the label	(E)	WEEE symbol (Explanation below)
	(B)	Invacare Logo	(F)	Type B Applied Part
	(C)	Serial Number	(G)	Recommendation to read the user manual before using the product.
	(D)	Reference Number		

This is the WEEE symbol (Waste Electrical and Electronic Equipment Directive).

This product has been supplied from an environmentally aware manufacturer. This product may contain substances that could be harmful to the environment if disposed of in places (landfills) that are not appropriate according to legislation.

- The ‘crossed out wheeled bin’ symbol is placed on this product to encourage you to recycle wherever possible.
- Please be environmentally responsible and recycle this product through your recycling facility at its end of life.

2 Setup



CAUTION!

Risk of Injury or Damage

If clamping lever or handle and housing screws are fitted too tightly or replaced by hardware not approved by Invacare, IDC internal parts might get damaged. This could lead to an unintended drive signal and the IDC could start running by itself when IDC is controlling the wheelchair.

- Do NOT tighten handle screws at all. Handle may get stuck or blocked.
- Tighten clamping lever and housing screws hand-tight only (max. 6 Nm).
- Tighten screw on monoblock to 1 Nm.
- Only use clamping lever and screws provided by Invacare for fixing IDC.
- Test power wheelchair before use.

2.1 Installing / Removing the IDC




NOTICE!

If an IDC is installed it will increase the footprint of the wheelchair. When driving without attendant, the IDC could get damaged, for example, when driving backwards against an obstacle.


- Remove or fold in the IDC before driving without attendant.

Installing




1. Position the adapter bracket (A) onto the KLICKfix adapter (B).
2. Push the IDC downwards until it locks automatically.

Removing



1. Press the red KLICKfix button (A) and hold it.
2. Lift the IDC upwards to remove it from the adapter.

Folding in the IDC



1. Remove the IDC from the adapter.
2. Press the red KLICKfix button (A) and hold it.
3. Position the lower edge of the adapter bracket (B) inside the lower fixation of the adapter (C).
4. Release the KLICKfix button.

Preparing IDC for Transport in a Vehicle




WARNING!

Risk of Injury and Damage due to Loose Parts during Transport

For transport IDC with LiNX must be completely removed from power wheelchair and stored safely.

- Disconnect and take off IDC for transport.

1. Disconnect bus cable.
2. Press and hold KLICKfix button .
3. Take off IDC.
4. During transport of power wheelchair store IDC safely, for example in glove compartment.

2.2 Connecting IDC to Controls System



The IDC takes over control of the wheelchair and may only be used by an attendant. It is imperative that the wheelchair occupant himself is not allowed to use the IDC.



Fitting should only be carried out by an Invacare provider. When fitting, ensure that the bus cable is routed and fixed correctly to avoid damage to the bus cable.

1. Power down wheelchair controls.
2. Connect one end of bus cable to IDC, other end to free bus socket on wheelchair.



In case there is no spare bus socket on controls system, you can order a 4-way bus socket from your Invacare provider.

3. The IDC is now connected to the controls.




For further information, refer to LiNX Service Manual.

3 Usage

3.1 Power up IDC

1. Connect IDC.
2. Power up power wheelchair.
3. If flash code appears, restart.
4. If flash code stays, check flash code, refer to *4.2 Fault Codes and Diagnosis Codes, page 10*.

 If necessary, make firmware update of power wheelchair, install latest drive program, read error via LiNX Access iOS tool or PC tool.


3.1.1 Switching Control Between Attendant and Occupant


The remote that powers up the power wheelchair is in control of the power wheelchair.

Pressing the power button on a remote not in control, returns control from attendant to occupant and vice versa, unless other settings have been made.

3.2 Controlling the Power Wheelchair with IDC


1. Press power button on IDC to power up power wheelchair. The display on the remote illuminates. The power wheelchair is ready to drive.
2. Set driving speed, refer to *3.2.1 How to Steer the Power Wheelchair, page 8*.
3. Pull trigger lever to set power wheelchair in motion.


 After a restart, the power wheelchair always drives forwards.




You can set whether the power wheelchair should travel forwards or backwards by pressing the switchover button for forwards or reverse gear **A**. If reverse gear has been selected, the arrow **B** illuminates.

3.2.1 How to Steer the Power Wheelchair

 Make yourself familiar with steering the power wheelchair via IDC, so that no hazards are posed for the occupant or others.




The electronic system contained in the handle **A** registers the movement made by the attendant's hand. Slight pressure on the handle to the left or right is translated into a drive command for the drive wheels so that bends can be travelled without the need for much force.



Travel Direction


You steer the power wheelchair by pushing the handle to the sides. Whether you are travelling forwards or backwards, you must push the left-hand side of the handle to make the power wheelchair turn to the left. If you push the right-hand side of the handle, the power wheelchair will turn to the right.

The maximum speed of the power wheelchair can be set by the user to match the personal requirements and the surrounding conditions. The current maximum speed set is displayed on the power wheelchair speed display and can also be preset there. You can carry out further fine adjustments to the maximum speed on the IDC using the speed potentiometer.



Setting Maximum Speed

1. Set maximum speed with speed potentiometer (A).
2. Pull trigger lever to start driving.



Controlling Speed

The more the trigger lever is pulled, the faster the power wheelchair drives.




To brake quickly, simply let go of the handle. The power wheelchair brakes. If you pull the trigger lever with full force, the power wheelchair carries out an "Emergency Stop".

3.3 "Emergency Stop" Function




When using Emergency stop function while driving, the power wheelchair brakes harder, than when you just let go of the trigger lever.

There are three ways to emergency stop the power wheelchair.



If the power wheelchair needs to be stopped immediately due to a hazard situation, pull the trigger lever (A) firmly as far as possible. This causes an immediate forced drive switch-off. Once you release the trigger lever (A) and the attendant LED stops flashing, you can start to drive the power wheelchair again immediately.



Alternatively, press "forwards or reverse" push button to stop the power wheelchair. The attendant LED stops flashing and you can start to drive the power wheelchair again immediately.



Or, press the power button. The power wheelchair needs to power up, before you can drive the power wheelchair again.

4 Troubleshooting

4.1 Fault Diagnosis

If the electronic system shows a fault, use the following fault-finding guide to locate the fault.



Ensure that the drive electronics system is powered up before starting any diagnosis.


If the Status Display is OFF:

- Check whether the drive electronics system is powered up.
- Check whether all cables are correctly connected.
- Ensure that the batteries are not discharged.

If a Flash Code is Displayed by the Status Indicator:

- Proceed to the next section.

4.2 Fault Codes and Diagnosis Codes




The status indicator is yellow when remote is active. If there is a fault detected by the LiNX system, the status indicator Ⓐ flashes yellow.

The number of flashes indicates the type of fault.

The table below describes the fault indication, and a few possible actions that can be taken to rectify the problem. The actions listed are not in any particular order and are suggestions only. The intention is that one of the suggestions may help you clear the problem. If in doubt, contact your provider.

Flash Code	Fault Description	Possible Action
1	Remote fault	<ul style="list-style-type: none"> • Check cables and connectors. • Check other remotes, if installed. • Contact your provider.
2	Network or configuration fault	<ul style="list-style-type: none"> • Restart power wheelchair. • Check cables and connectors. • Recharge the batteries. • Check charger. • Contact your provider.
3	Motor 1 ¹ fault	<ul style="list-style-type: none"> • Check cables and connectors. • Contact your provider.
4	Motor 2 ¹ fault	<ul style="list-style-type: none"> • Check cables and connectors. • Contact your provider.
5	Motor 1 ¹ magnetic brake fault	<ul style="list-style-type: none"> • Check cables and connectors. • Check left magnetic brake is engaged. • Refer to the chapter “Pushing the Power Wheelchair in Freewheel Mode” in the user manual of your power wheelchair. • Contact your provider.
6	Motor 2 ¹ magnetic brake fault	<ul style="list-style-type: none"> • Check cables and connectors. • Check right magnetic brake is engaged. • Refer to the chapter “Pushing the Power Wheelchair in Freewheel Mode” in the user manual of your power wheelchair. • Contact your provider.
7	Module fault (other than remote module)	<ul style="list-style-type: none"> • Check cables and connectors. • Check modules. • Recharge batteries. • If the power wheelchair was stalled, reverse away or remove obstacle. • Contact your provider.

1 Configuration of the motors depending on the power wheelchair model



The attendant in control LED is green when IDC is in charge. If there is a fault detected by the LiNX system, the attendant ② in the control LED flashes green.

Flash Code	Fault Description	Possible Action
Flashing constantly	OON	Release trigger lever. If fault code stays: <ul style="list-style-type: none"> • Do not use IDC. Replace IDC. • Contact your provider.

5 Maintenance

General information about maintaining and cleaning remote and power wheelchair are described in the user manual of the power wheelchair.

Inspection Checks

The following table lists inspection checks that should be performed by the user and their intervals. If the power wheelchair fails to pass one of the inspection checks, refer to the chapter indicated or contact your authorised Invacare provider. A more comprehensive list of inspection checks and instructions for maintenance work can be found in the service manual for the corresponding wheelchair, which can be obtained from Invacare. That manual, however, is intended to be used by trained and authorised service technicians, and describes tasks which are not intended to be performed by the user.

Item	Inspection Check Before Each Use with IDC	If Inspection is not Passed
IDC	Check for damage and wear.	<ul style="list-style-type: none"> • Check all functionalities below.
IDC	Check if trigger lever for driving / emergency stop is completely released.	<ul style="list-style-type: none"> • Immediately stop using the IDC. • Contact your provider.
IDC	Check emergency stop functions.	<ul style="list-style-type: none"> • Immediately stop using the IDC. • Contact your provider.
IDC	Check power button for correct functionality.	<ul style="list-style-type: none"> • Contact your provider.
IDC	Check of correct operating mode for all LED flash lights.	<ul style="list-style-type: none"> • Contact your provider.
Electronics / connectors	Check all cables for damage and all connecting plugs for snug fit.	<ul style="list-style-type: none"> • Contact your provider.
Housing	Check IDC housing for loose parts.	<p>If there is a rattling noise, when you shake it:</p> <ul style="list-style-type: none"> • Immediately stop using the IDC. • Replace IDC. • Contact your provider.

**Australia:**

Invacare Australia Pty. Ltd.
Unit 18/12 Stanton Road,
Seven Hills, NSW 2147,
Australia
Phone: 1800 460 460
Fax: 1800 814 367
orders@invacare.com.au
www.invacare.com.au

EU Export:

Invacare Poirier SAS
Route de St Roch
F-37230 Fondettes
Phone: (33) (0) 2 47 62 69 80
serviceclient_export@invacare.com
www.invacare.eu.com

Canada:

Invacare Canada LP
c/o Motion Concepts Canada
84 Citation Drive Unit 11
Concord, Ontario L4K 3C1

New Zealand:


Invacare New Zealand Ltd
4 Westfield Place, Mt Wellington
1060
New Zealand
Phone: 0800 468 222
Fax: 0800 807 788
sales@invacare.co.nz
www.invacare.co.nz

United Kingdom & Ireland:

Invacare Limited
Pencoed Technology Park, Pencoed
Bridgend CF35 5AQ
Phone: (44) (0) 1656 776 200
uk@invacare.com
www.invacare.co.uk

UKRP Invacare UK Operations Limited
Unit 4, Pencoed Technology Park, Pencoed
Bridgend CF35 5AQ
UK

1643411-C 2024-10-10



Making Life's Experiences Possible®



Yes, you can.®