



INVACARE® TOP END®
FORCE™ K
HANDCYCLE



The **Invacare® Top End® Force™ K Handcycle** is a kneeling position handcycle available in two custom, made-to-measure aerodynamic designs. The Top End Force K Kneeler is used by the most able-bodied athletes, while the Top End Force K Double Amputee Kneeler is custom designed for the double above-the-knee amputee. What makes the Top End Force K Handcycle different from the rest of the Invacare Top End Force Series Handcycles? The Top End Force K Handcycle is a kneeling-only, non-adjustable position handcycle that requires dialed-in measurements and engineering drawings - all of which are easy to find on the website: www.topendwheelchair.com. Can you force the pace?


Features

7005 aluminium custom welded frame
7005 aluminium aerodynamic fork
Breathable mesh upholstery
Seat and back cushions
Tapered crank arms
Frictionless pedal bearings
Aluminium ergonomic hand pedals
SRAM® 30 speed drive train
Chainrings: 53-39 or 48-39-24
SRAM cassette: 11-34

Tektro® Caliper Brake
Fork mounted locking second brake
Steering dampener
26" high performance wheels
3°, 6°, 9° or 11° of camber
Threaded axles
Drafting bumper
Parking brake
Safety flag
Chain guard

Many paint choices
Reflective stickers
250 lb. weight capacity

Model #: FRCK



Invacare Corporation
www.invacare.com

USA

One Invacare Way
Elyria, Ohio
44035-4190
(800) 333-6900

Canada

570 Matheson Blvd. E.,
Unit 8
Mississauga, Ontario
L4Z 4G4 Canada
(800) 668-5324

Top End

www.topendwheelchair.com
Sports and Recreation Products
4501 63rd Circle North
Pinellas Park, Florida
33781-5914
(727) 522-8677
(800) 532-8677
Fax: (727) 522-1007

©2014

Invacare Corporation.
All rights reserved.
Trademarks are identified
by the symbols™ and
®. All trademarks are
owned by or licensed to
Invacare Corporation unless
otherwise noted. SRM is a
trademark of SRM Service
Center Inc. USA. SRAM is a
registered trademark of SRAM
Corp. Tektro is a registered
trademark of Tektro USA.
Specifications are subject to
change without notification.
Form No. 12-451
Rev. 12/15
141159

Options

Carbon fiber wheels
Component, brake and cassette upgrades
SRM® power meter measuring tool integrated
directly into Top End cranks

See the Force K Handcycle order form for all of the
great options available!

Frame Styles:

Invacare Top End Force K Kneeler
designed for the most able athletes with a custom
knee/seat bucket.

Invacare Top End Force K Double Amputee Kneeler
designed for the double above-the-knee amputee
utilising a custom cage.

