

Comet[®] Ultra



Built to support with performance
and reliability



Ideal for those with larger body shapes, the **Comet^{Ultra}** offers a secure, supportive and powerful driving experience without compromising on style or functionality. Packed with all of the safety features and driving characteristics seen throughout the **Comet** family, the **Comet^{Ultra}** will get you to your destination quickly and efficiently, in complete comfort and control.

Features and options



Improved suspension

Improved suspension to offer superior comfort to larger individuals.



Reinforced seat column

Offers additional stability and safety when seated



Reinforced height and width adjustable armrest

Can be easily adjusted to the user and offers additional support when getting on and off the scooter

Power, reliability and customisation









Due to a carefully thought out design process and a robust construction, the safety and stability of the **Comet^{Ultra}** is never a concern; even for users weighing up to 220kg. Developed with the highest quality components and safety features, the **Comet^{Ultra}** makes a trustworthy companion no matter where you go. Furthermore, its flexibility and adaptability means you can customise your scooter to suit your unique style and individual preference.









Complete comfort and support


The special contoured Ultra Seat has been developed in conjunction with therapists to ensure an optimum seating position and increased comfort. In addition, a pneumatically adjustable lumbar support allows for further comfort when seated if and when required. The reinforced seat base and adjustable armrests offer excellent support and allows the use of the armrest as additional support when getting on and off the scooter.

Technical data

For more comprehensive pre-sales information about this product, including the product's user manual, please see your local Invacare website.

	Backrest width	Seat depth	Seat to floor height	Backrest height	Driving unit width	Total length	Empty weight	The heaviest part
Comet^{Ultra}	 660 mm	 535 mm	 Min. 440 mm Max. 510 mm	 515 mm	 665 - 685 mm	 1450 mm	 148 kg	 approx. 63 kg

	Max. user weight	Battery capacity	Turning diameter	Max. climbable obstacle height	Ground clearance	Max. safe slope	Max. driving distance	Max. driving speed
Comet^{Ultra}	 220 kg	 2 x 12 V/75 Ah (C20) AGM	 2750 mm	 100 mm	 100 mm	 1°	 49 km	 10 km/h

	Motor capacity
Comet^{Ultra}	 10km/h: 650 W / 1800 W ¹

Colour of shrouds

Marble White	Alium Blue	Jade Green	Luxor Beige	Crystal Copper
				
Carmine Red	Iconic Chrome	Diamond Silver	Electric Blue	Garnet red
				

Upholstery

Black Vinyl


1. Peak performance

Tel: +41 61 487 70 70
hqeurope@invacare.com
www.invacare.eu.com

Invacare International GmbH
Benkenstrasse 260
4108 Witterswil
Switzerland

