



Küschall® CHAMPION SL

MADE TO
TRAVEL

—
BUILT TO
PERFORM



Kuschall.com



*FOLDABLE
AND COMPACT*

CHAMPION

SL

Meet the all new Kuschall Champion SL, an exciting and stylish new Kuschall design that retains all the travel-conscious characteristics of the original Champion series but welcomes a lighter and more agile performance. Engineered with the increased stiffness and rigidity of 'The KSL', the new Champion SL gives optimum driver control and high levels of responsiveness for an even greater drive.

*SUPER LIGHT
AND RIGID*





The Champion SL promises to give the most discerning of active wheelchair users an impeccable drive.



We poured our passion and expertise into every detail of the Champion SL.

Our goal was to design a compact, foldable chair that gives an exceptional drive. Not only did we commit to the highest standards of functionality and durability in our design, but we reached for a chair that captures the spirit and aesthetic of Kuschall.



MADE TO TRAVEL *BUILT TO PERFORM*

Designed for an active lifestyle, the Champion SL helps take you on a journey to explore new adventures, pursue passions and help unlimit your world.





TRANSPORTS WITH EASE

We would never change what you love about the Kuschall Champion. A travel-conscious design, the Champion SL is a winner when it comes to convenience. Simply fold the frame for easy transportation and compact storage. The new Champion SL takes you places you want to go – with confidence.





*EXCEPTIONAL
DRIVE
PERFORMANCE*

The Champion SL is a highly versatile active chair built especially to perform. Super light and robust, we have focused on a design precisely engineered for energy efficiency. With increased resistance and responsiveness, the Champion SL pays attention to efficient propulsion and optimal biomechanics.





TAKE A CLOSER LOOK

SUPER LIGHT

Now 5.9kg transport weight for easier handling.

RESPONSIVE

Gives easy and quick changes in direction, speed and navigation of tight corners.

TRANSPORTS WITH EASE

Unique folding mechanism makes transporting and storage simple.



ERGONOMIC GRIP

For easier repositioning and transfers.

HYDROFORMED TECHNOLOGY

Strengthens the frame for increased rigidity and improved energy efficiency when propelling.



**WATCH
THE
VIDEO**



CUSTOMISE YOUR CHAMPION SL

1. FRAME COLOUR

Standard frame colours



Pure White Anthracite Grey Space Black



Carmine Red Azure Blue

Extra frame colours*



Flashy Yellow Black Glossy Silver Grey Royal Blue



Ruby Red Dynamite Orange Forest Green

*Extra lead time


2. DECAL STYLE



3. DECAL COLOUR



Red Black



Light Grey Gold Red Blue Dark Grey Black

4. UNIQUE SERVICE

Create your new Küschall wheelchair from the wheels up with UNIQUE. This dedicated service will create your personalized chair, designed to fit you and your lifestyle.

UNIQUE
Designed by you. Built by us.



Pictured: Spinergy wheel with UNIQUE features. Scan the QR code on this page for more information.





TECHNICAL SPECIFICATION

SEAT WIDTH

360 - 480 mm
In increments of 20 mm



SEAT DEPTH

340 - 480 mm
In increments of 20 mm



KNEE-TO-HEEL LENGTH

320 - 500 mm
Stepless adjustable



BACKREST HEIGHT

300 - 465 mm
In increments of 15 mm



BACKREST ANGLE

76° / 80.5° / 85° / 89.5° / 94°



SEAT HEIGHT

Front: 450 - 540 mm
Rear: 390 - 490 mm



FRAME ANGLE

75° / 90°
Straight or tapered



CENTER OF GRAVITY

8 positions



CAMBER

0° / 1° / 3°



TOTAL LENGTH

75° approx. 930 mm
90° approx. 830 mm



TOTAL WIDTH

SW + 160 - 240 mm



TOTAL WIDTH FOLDED

Approx. 280 - 340 mm



TOTAL WEIGHT UNLOADED

Approx. 8.4 kg



TRANSPORT WEIGHT

Approx. 5.9 kg



MAX. USER WEIGHT

100 kg



All measurements and weights indicated are based on seat width 400 mm and seat depth 400 mm in the lightest and most compact configuration. Weight and measurements can change with the different configurations. For further information please check the user manual on our website.



UNLIMIT YOUR WORLD

Tel: +41 61 487 70 70
info@kuschall.com
www.kuschall.com

Invacare International GmbH
Neuhofweg 51
4147 Aesch BL
Switzerland

