

PRECISION • COMFORT • SUPPORT



Tel: +41614 877070  
hqeurope@invacare.com  
www.invacare.eu.com



**matrx<sup>®</sup>MAC**



NEW

A multi-adjustable contour back designed for more complex postural needs in tilt and recline applications. It features three lateral depths that are interchangeable and ambidextrous to mix and match for bespoke trunk support.



**matrx<sup>®</sup>HUG**



NEW

Optional cushion insert with Stamina Fibre<sup>®</sup> for maximum immersion and envelopment. Microclimate control, hypoallergenic and machine washable. Available on Matrx MAC back and Matrx headrest pads.



**matrx<sup>®</sup>FLIPP**



NEW

Optional flip back feature for easy transfers, that doesn't affect the initial set up of the headrest. Available on Matrx Elan and Matrx Loxx headrest.



PRECISION • COMFORT • SUPPORT

We are committed to enhancing people's lives through our clinician-led seating solutions.



**matrx<sup>®</sup>**  
SCIENCE OF BEATING

THE RANGE  
2025





## matrx<sup>®</sup>LIBRA



Our signature modular positioning cushion that focuses on immersion, envelopment and offloading. It provides our highest level of skin protection and stability.



## matrx<sup>®</sup>MX2



The carefully crafted carbon fibre design makes this super-lightweight and shock-absorbing back, a solution for active wheelchair use.



## matrx<sup>®</sup>E2



Our signature ergonomic back available in 2 depths for optimal support and pressure distribution. Featuring Ripple Foam technology for microclimate control and enhanced immersion.



## matrx<sup>®</sup>PB



Contoured back providing exceptional midline support. Featuring durable moulded foam which wraps around the backshell for optimal skin protection.



## matrx<sup>®</sup>LOXX



A simple but highly adjustable headrest with precise positioning for tilt-in space wheelchairs. With intuitive, four points of adjustment and offsetting capability for a comfortable fit.



## matrx<sup>®</sup>ELAN



Award winning, highly adjustable headrest for more complex positioning. The single tool hardware offers simple set up. Six points of adjustment and offsetting capability allows for the perfect fit.