

AQUATEC®

Bath



Folding shower seat &
bathboard



Yes, you can.®



Aquatec® Sansibar

- Space-saving: mounted on the wall and can be folded when not in use
- Large hygiene recess makes personal hygiene easy
- Height-adjustable support legs for individual adjustment and additional stability
- Compensates for floor unevenness through finely adjustable feet

Invacare Marina


- Attractively designed bathtub board
- Non-slip seat surface
- Hand grip provides additional safety
- Practical shower head holder on both sides
- Individually adjustable using adjustment elements
- Also available as an XL version

Technical data


For more comprehensive pre-sales information about this product, including the product's user manual, please see your local Invacare website.

Order number	Aquatec® Sansibar	Invacare Marina	Invacare Marina XL
Folding shower seat	9.10.101		
Bathtub board with hand grip		1541330	
Bathtub board XL			1471594
Dimensions			
Seat width	400 mm (A)	-	-
Seat depth	410 mm (B)	-	-
Seat height	460 - 530*1 mm (C)	-	-
Total width	400 mm (D)	360 mm (A)	350 mm (A)
Total depth	460 mm (E)	690 mm (B)	740 mm (B)
Total height	530 - 600 mm (F)	160 mm (C)	160 mm (C)
Width, hand grip	-	180 mm (D)	180 mm (D)
Adjustment range, inner dimension	-	400 - 620 mm (E)	420 - 690 mm (E)
Graduation, height adjustment	23 mm	-	-
Depth, folded	110 mm	-	-
Height, folded	540 - 610 mm	-	-
More info			
Weight	2.6 kg	1.6 kg	2.8 kg
Load capacity	130 kg	150 kg	200 kg
Colour	White	White	White

Aquatec Sansibar



Invacare Marina / Marina XL



1. 3. Height adapter optional

Tel: +41 61 487 70 70
 hqeuropa@invacare.com
 www.invacare.eu.com

Invacare International GmbH
 Neuhofweg 51
 4147 Aesch BL
 Switzerland



Yes, you can.®